



OMEGA LA1 WHEELCHAIR

Features:

- lightweight aluminium frame
- swing away removable foot rests
- fold down back rest
- 8" rear wheels, attendant brakes
- easy to fold for travel or storage
- user weight capacity of 125kg



OMEGA LA1 WHEELCHAIR	Metric	Imperial
Approximate Total Product / Unit Weight	8.9 kg	19.6 lbs
Overall Width (product or frame when opened)	640 mm	25.2 inch
Overall Depth (product or frame when opened)	1,060 mm	41.7 inch
Width / Depth (product or frame when closed)	290 mm	11.4 inch
Overall Height (or height range if adjustable)	940 mm	37.0 inch
Armrest Width (outside dimension)	560 mm	22.0 inch
Armrest Width (inside dimension)	440 mm	17.3 inch
Armrest Height (from top of seat)	220 mm	8.7 inch
Back / Back-Rest Height (from top of seat)	410 mm	16.1 inch
Seat Height (from top of seat to floor)	470 mm	18.5 inch
Seat Width (actual seat, not available space)	460 mm	18.1 inch
Seat Depth (actual seat, not available space)	400 mm	15.7 inch
Wheels (dimensions of front)	152 mm	6.0 inch
Wheels (dimensions of rear)	203 mm	8.0 inch
Recommended Maximum User Weight	125.0 kg	275.6 lbs



Max Mobility
..... live life to the max