

## EZ FOLD-N-GO ROLLATOR

### Features:

- portable - quickly & easily folds up for storage in small spaces
- 6 inch allows for trouble-free maneuvering over all surfaces
- convenient seat - perfect for resting
- height adjustable - and supports up to 115kg
- contemporary look - three stylish colours to choose from



Black Walnut



Regal Rose



Cobalt Blue



FOLD-N-GO WALKER	Metric		Imperial	
Collapsed Dimensions	90.8 x 33.78 x 25.4	cm	35.75x13x10	nch
Wheels	15.24	cm	6	nch
Seat Dimension	46.99 x 12.7	cm	18.5W x 5 D	inch
Height Adjustment	83.82 x 91.44	cm	33 x 36	inch
User Height	147.32 to 203.2	cm	4'10" to 6'8"	
Open Depth	58.42	cm	23	inch
Open Width	69.85	cm	27.5	inch
Open Handle Width	13.18	cm	17	inch
Weight	6.35	kg	14	lbs
Weight Capacity	115	kg	250	lbs