


FOLD-N-GO WALKER

Features:

- Portable - quickly & easily folds up for storage in small spaces
- 6 inch swivel wheels with locking pin and rear Easy-Glide feet - allows for trouble-free maneuvering
- height adjustable - and supports up to 140kg
- contemporary look - three stylish colours to choose from



Black Walnut



Cobalt Blue



Regal Rose



FOLD-N-GO WALKER	Metric		Imperial	
Collapsed Dimensions	91.44 x 13.3 x 20.32	cm	36 x 5.25 x 8	inch
Wheels	15.24	cm	6	inch
Height Adjustment	81.28 x 97.79	cm	32 x 38.5	inch
User Height	147.32 to 203.2	cm	4'10" to 6'8"	
Open Depth	54.61		21.5	inch
Open Width	59.69		23.5	inch
Open Handle Width	45.72		18	inch
Weight	3.4	kg	7.5	lbs
Weight Capacity	140	kg	287	lbs

