





ZETA PE1 PEDAL EXERCISER

Features:

- can provide gentle exercise
- for use on the floor or a table
- compact & easy to fold for travel
- tension adjustable to suit user
- flexible feet to minimise slipping
- unique multi-function counter



ZETA PE1 PEDAL EXERCISER	Metric	Imperial
Approximate Total Product / Unit Weight	2.0 kg	4.4 lbs
Overall Width (product or frame when opened)	385 mm	15.2 inch
Overall Depth (product or frame when opened)	460 mm	18.1 inch
Width / Depth (product or frame when closed)	270 mm	10.6 inch
Overall Height (or height range if adjustable)	250 mm	9.8 inch