

CRYSTAL² SELF PROPEL WHEELCHAIR

RG18B

**NEW
PRODUCT**

FEATURES

Lightweight: 12.5kg (Aluminium frame)
Compact - Folding frame and fold down backrest
Removable / height adjustable swing away legrests
Fold up footplates with heel straps
Attendant operated hand / park brakes
Occupant operated park brakes
Padded armrests
Quick release self propel rear wheels

Quick release wheels

Swing away legrests

SPECS

(Dims shown are nominal)

Overall (Including swing away legrests):

660mm (W) x 1080mm (L) x 925mm (H)

Folded Overall (Legrests / rear wheels removed):

300mm (W) x 670mm (L) x 735mm (H)

Seat: **450mm (W) x 420mm (D) x 520mm (H)**

Armrest: **300mm (L) x 50mm (W) x 220mm (H)**

Backrest: **450mm (W) x 400mm (H)**

Footrest: **175mm (W) x 155mm (D) x 340 - 440mm (H)**

Wheels: **Ø200mm / 8" (Front) - Ø610mm / 24" (Rear)**

Product Weight: **12.5kg**

Max. User Weight: **127kg**

IMPORTANT NOTE:

It is the responsibility of the user to ensure the suitability of this product for its intended use.

The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.

Please contact **AbleMedilink** if you require more information on this product.

P: 03 9557 9038

5 Nicholson St, Bentleigh VIC 3204, Australia

ablemedilink.com.au