

“CrossOver” Seat Walker & Occasional Transit Chair Walker

RG4408



Can be used with
or without
removable
leg rests



Features

- Side folding seat walker
- Converts to a transit wheelchair
- Aluminium construction
- Top quality powder coated finish
- Side folding action
- Soft operation hand brake
- Seat height at 560mm from floor
- 115mm wide padded backrest
- Removable swing away leg rests fitted
- Walking stock holder fitted as standard
- 200mm spoked wheels
- Available in three colours
Blue / Red / White

Dimensions

Seat Height	560mm
Seat Width	450mm
Seat Depth	310mm
Unit Overall Width	590mm
Unit Overall Length	980mm
Handle Height Range	900mm to 1000mm
Unit Net Weight	9.42kg
Unit Net Weight no footrests	8.3kg
Wheel Size	200mm / 8"
Maximum User Capacity	120kg

Phone: 03 9557 9038
5 Nicholson St, Bentleigh VIC 3204
ablemedilink.com.au



REDGUM Brand is a 100% Australian Owned & Operated