

Tri Walker with Basket ,Tray & Zip Bag



Features

- Strong steel construction
- Powder coated folding frame
- Multi height adjustable handles
- 360degree swivelling front wheels
- 195mm Micro cellular castors
- Non marking castors
- Locking handbrakes
- Fitted with soft black zip bag
- Fitted with powder coated basket
- Fitted with removable tray
- Available in either Blue or Red

Dimensions

- Maximum user capacity 100kg
- Maximum width 660mm
- Overall length 580mm
- Height adjustment 880 to 980mm
- Net weight 6.5kg (w/o basket)

Phone: 03 9557 9038
5 Nicholson St, Bentleigh VIC 3204
ablemedilink.com.au



REDGUM Brand is a 100% Australian Owned & Operated