

## Pacer Aluminium Seat Walker



**PACER 6" Seat Walker**  
Model RG4205AB- Hand brake 56.5cm  
Others Available:  
Model RG4202AC- Push down brake  
Model RG4208AJ – Hand brake 48cm seat height

### Dimensions

- Seat height 540mm
- Max user capacity 120kg
- Walker net weight 5.5kg – 7.5kg



**PACER 8" Seat Walker**  
Model RG4206AC – Hand Brake 56cm seat height  
Other Available:  
Model RG4206AB - Hand Brake 51cm seat height  
Model RG4206AD - Hand Brake 61cm seat height

### Features

- Aluminium powder coated frame
- Available in three marble finish colours
- Multi height adjustment handles
- Large soft under seat storage bag
- Comfortable folding backrest

---

Phone: 03 9557 9038  
5 Nicholson St, Bentleigh VIC 3204  
[ablemedilink.com.au](http://ablemedilink.com.au)



*REDGUM Brand is a 100% Australian Owned & Operated*