

NEW  
PRODUCT


# BLANKET SUPPORT

RG64BS

SUPPORTS BEDDING  
above person in bed



FOLDS FLAT  
for storage and transport



## FEATURES

- Supports bedding above person in bed
- Two width settings: 350mm / 460mm
- Easy to remove and fit when making bed.
- Folds flat for transport and storage (60mm)
- Lightweight (2kg)

## SPECIFICATIONS

Overall: **510mm L x 460mm W x 560mm H**

Overall: **560mm L x 460mm W x 60mm H**  
(Folded)

Product weight: **2kg**

Max. Supportive Load: **10kg**

Please contact **REDGUM Brand** if you require further information on this product

P: **08 9248 4180**

E: **orders@for-de.com.au**

A: **1 Business way Malaga Western Australia 6090**

Part of the **For-de Group Pty Ltd**

ASK FOR THE BRAND

redgumbrand.com.au