

NEW PRODUCT



**LIMITED
EDITION**
**METALLIC
SKY BLUE**



**FOLDS UP
EASILY FOR
TRANSIT OR
STORAGE**

IMPORTANT NOTE:

It is the responsibility of the user to ensure the suitability of this product for its intended use. The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.

EASY FOLD SEAT WALKER

RG4403

FEATURES

- Height adjustable handles
- Ergonomic rubber handgrips
- User friendly loop brake handles
- Brake + park brakes
- Mesh seating
- Padded Backrest
- Walking stick holder
- Removable storage bag
- Side and front reflectors
- Large wheels for excellent mobility
- Folds up easily for transit or storage

SPECIFICATIONS

(All dimensions shown are nominal)

- Overall: **700mm (L) x 680mm (W) x 880mm (H)** (at lowest setting)
- Overall folded: **780mm (L) x 230mm (W) x 280mm (H)**
- Handle Height: **880 - 980mm** (5 settings)
- Seat: **480mm (W) x 240mm (D) x 550mm (H)** (floor to seat)
- Wheel Diameter: **Ø200mm / 8"**
- Carry Bag: **370mm (W) x 220mm (D) x 200mm (H)**
- Product Weight: **8.5kg**
- Maximum User Weight: **136kg**
- Maximum Bag Loading Weight: **5kg**

Please contact **REDGUM Brand** if you require further information on this product

P: 08 9248 4180

E: orders@for-de.com.au

1 Business Way Malaga Western Australia 6090
Part of the **For-de Group Pty Ltd**

ASK FOR THE BRAND

redgumbrand.com.au